

NANOFIL[®]

A L W A Y S F O R Y O U R S M I L E

Dental Filling Resin

Light Curing Composite - Flowable

Scope

- Class III and V restorations
- Restoration of minimally invasive non stress-bearing cavity
- Base/liner under direct restorations
- Repair of small defects in esthetic indirect restorations
- Pit and fissure sealant
- Undercut blockout
- Repair of resin and acrylic temporary materials

Specifications

3g/syringe

Shades

A1, A2, A3, A3.5, B1, B2



ATBIO
AT&M Biomaterials Co., Ltd.



☎ 010-56764321

🌐 www.atmbio.com

✉ info@atmbio.com

📍 No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.

NANOFIL[®]

A L W A Y S F O R Y O U R S M I L E

VERSATILITY OF USE

High-quality products new energy, a wide range of use

VARIOUS COLORS AVAILABLE

Perfectly match tooth color

LOW POLYMERISATION SHRINKAGE


Polymerisation shrinkage rate <2.4%

EXCELLENT FLOWABILITY

Full penetration, extension effect, adapts to small crevices.

STRONG ANTI-WEAR ABILITY

Closer to the original condition of the teeth.



CE 0197

ATBIO
AT&M Biomaterials Co., Ltd.



010-56764321



www.atmbio.com



info@atmbio.com



No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.