

NANOFIL[®]

A L W A Y S F O R Y O U R S M I L E

Dental Pit&Fissure Sealant

Scope


Pit & Fissure Sealant is designed for sealing the enamel pits and fissures of teeth to aid in the prevention of carries.

Specifications

3g/syringe

Color

White 、 blue



☎ 010-56764321

🌐 www.atmbio.com

✉ info@atmbio.com

📍

No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.

CE₀₁₉₇

ATBIO

AT&M Biomaterials Co., Ltd.

NANOFIL[®]

A L W A Y S F O R Y O U R S M I L E

RAPID CURE

The curing time is only 20 seconds

LOW VISCOSITY

Viscosity < 2000cp

CONVINIENT OPERATION

Needle design, easy coating

TWO-TONE DESIGN

White and blue are available

GOOD STRENGTH

Flexural strength ≥ 130 and
compressive strength ≥ 120 after curing

CLINICAL VERIFLCATION

Has the excellent pit and fissure sealing effect



010-56764321



www.atmbio.com



info@atmbio.com



No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry
Base, Haidian District, 100094 Beijing, P.R. China.



ATBIO

AT&M Biomaterials Co., Ltd.